

# TECHNICAL DATASHEET

> Protective Equipment

## Vogt Medical Disinfectant



**Intended Use:** Hand and surface disinfectant. Disinfectant recommended by the World Health Organization. To combat the spread of COVID-19.

### Product Description

- > For consumer and professional use.
- > The product contains ethanol, which provides its antibacterial and viricidal properties.
- > Application for hand disinfection: apply a small amount (3-5 ml) of the product to the target area and rub it well over your skin, leave for at least 15-30 seconds until it dries.
- > Application for surface disinfection: for rapid disinfection in the medical industry, nursing homes, sanatoriums, fitness and SPA centres, the beauty industry, kindergartens and schools.
- > Use the spray method to disinfect surfaces. Pour into the sprayer, spray the surface and wait for it to dry. If surfaces are dirty, they must be cleaned before disinfection.
- > Frequency of use as necessary.

---

### TD > Protective Equipment

created on: June 2020 by: Liudmyla Honchar

approved on: July 2020 by: Liudmyla Honchar

# TECHNICAL DATASHEET

## VOGT MEDICAL DISINFECTANT

### PRODUCT FORMULATION

Type of formulation: Liquid

Active substance: Ethanol 80%

Ingredients: Ethanol 96%, Aqua, Hydrogen Peroxide, Glycerin

- Ethanol provides anti-bacterial and viricidal properties.
  - Glycerol: used as a moisturizer, but other softeners may also be used for skin care, provided that they are soluble in water and alcohol and do not increase toxicity or promote allergies.
  - Hydrogen peroxide: used for inactivation of microbial spores that contaminate the solution, and it is not an active substance for hand disinfection.
- 

### QUALITY AND SAFETY

Conforms to Regulation (EC) No. 1907/2006 (REACH), Annex II

#### STORAGE

Temperature

Between +5°C +25°C

#### PACKAGING

Unitary

Volume: 100 ml, 150 ml, 200 ml, 250 ml, 500 ml, 1l, 5l

---

#### TD > Protective Equipment

created on: June 2020 by: Liudmyla Honchar

approved on: July 2020 by: Liudmyla Honchar