

UROLOGY PRODUCTS



DE VM® Foley Katheter

- › Zeichnet sich durch einen hohen Tragekomfort, individuelle Größen und hochwertige, innovative Materialien
- › Es ist einfach von 3- auf 2-lumigen Katheter zu wechseln, ohne den Anbieter wechseln zu müssen, das bedeutet, dass sich der Patient nicht an einen anderen Katheter eines Mitbewerbers gewöhnen muss und ihm so der gewohnte Tragekomfort nicht verloren geht
- › Werden zur dauerhaften Harnableitung eingesetzt
- › Ein transparenter und symmetrischer Ballon garantiert maximalen Schutz bei der Verwendung
- › Silikonbeschichtung verhindert Reizung der Harnröhre

Materialien:

Latex mit Silikonbeschichtung / 100% Silikon

Sterilisation:

Ethylenoxid

Weg:

2 oder 3

Größe (CH/FR):

- für Kinder: 8 – 10
- für Erwachsene: 12 - 26

Länge (mm):

280; 400

Max. Volumen-Ballon (ml):

3 – 30

Produktnormen:

EN 1616

Risikoklasse:

Ila / IIb

Verpackung:

Blister / Box / Karton

Haltbarkeit:

5 Jahre

EN VM® Catheter Foley

- › Are designed to be used for a long-term urinary drainage
- › Transparent, symmetrical balloon guarantees maximum protection
- › Two large, lateral eyes ensure optimal urinary flow
- › The product range includes 2-way or 3-way catheters
- › Latex catheters bring maximum comfort to patients
- › Silicone coating prevents irritation of the urethra and reduces the risk of an allergic reaction
- › The colour code on the connector and the fill level volume indicator for the catheter balloon guarantee an additional convenience

Materials:

Latex with silicone coating / 100% Silicone

Sterilization:

Ethylene oxide

Way:

2 or 3

Size (CH/FR):

- for children: 8 – 10
- for adults: 12 – 26

Length (mm):

280; 400

Max. volume-ballon (ml):

3 – 30

Compliance standards:

EN 1616

Risk classification:

Ila / IIb

Packaging:

blister / box / carton

Shelf life:

5 years

VM® Catheter Foley

Index No.	Material	Way	Size (CH/FR)	Length (mm)	Max. volume balloon (ml)	Colour
894	latex silicone coated*	2	8	280	3-5	black
895	latex silicone coated*	2	10	280	3-5	grey
896	latex silicone coated	2	12	400	5-10	white
897	latex silicone coated	2	14	400	5-10	green
898	latex silicone coated	2	16	400	5-10	orange
899	latex silicone coated	2	18	400	5-10	red
900	latex silicone coated**	2	20	400	5-15	yellow
901	latex silicone coated**	2	22	400	5-15	violet
902	latex silicone coated	2	12	400	30	white
903	latex silicone coated	2	14	400	30	green
904	latex silicone coated	2	16	400	30	orange
905	latex silicone coated	2	18	400	30	red
906	latex silicone coated	2	20	400	30	yellow
907	latex silicone coated	2	22	400	30	violet
908	latex silicone coated	2	24	400	30	blue
909	latex silicone coated	2	26	400	30	pink
910	100% silicone*	2	8	280	3-5	black
911	100% silicone*	2	10	280	3-5	grey
912	100% silicone	2	12	400	5-10	white
913	100% silicone	2	14	400	5-10	green
914	100% silicone	2	16	400	5-10	orange
915	100% silicone	2	18	400	5-10	red
916	100% silicone	2	20	400	5-10	yellow
917	100% silicone	2	20	400	5-15	yellow
918	100% silicone	2	22	400	5-10	violet
919	100% silicone	2	22	400	5-15	violet
920	100% silicone	2	24	400	5-10	blue
921	100% silicone**	2	18	400	20	red
922	100% silicone	2	20	400	30	yellow
923	100% silicone	2	22	400	30	violet
924	100% silicone	2	24	400	30	blue
925	latex silicone coated	3	16	400	30	orange
926	latex silicone coated	3	18	400	30	red
927	latex silicone coated	3	20	400	30	yellow
928	latex silicone coated	3	22	400	30	violet
929	latex silicone coated	3	24	400	30	blue
930	latex silicone coated	3	26	400	30	pink
931	100% silicone	3	16	400	5-10	orange
932	100% silicone	3	18	400	5-10	red
933	100% silicone**	3	18	400	20	red
934	100% silicone	3	20	400	30	yellow
935	100% silicone	3	22	400	30	violet
936	100% silicone	3	24	400	30	blue

*Pediatric

**not for sale in RU