

GYNECOLOGY PRODUCTS



DE VM® Gynäkologische Kits

- › Wurde für gynäkologische Untersuchungen und geburts- hilffiche und gynäkologische Verfahren entwickelt
- › Sind in verschiedenen Konfi- gurationen erhältlich und kön- nen durch die Kombination der einzelnen Produkte nach den individuellen Vorlieben angepasst werden:
 - Vaginal Spekulum
 - Schuhüberzüge
 - Serviette
 - Gynäkologischen Kollektoren
 - Untersuchungshandschuhe
 - Kondom
 - Mikroskop-Objekträger

Materialien:

PS; PE; PP; Vliesstoff; Latex / Nitril; Nylon; Stahl; ABS

Sterilisation:

Ethylenoxid

Kit Typ:

MINI; BASIC; CYTO-Standart; CYTO-Premium; ENDO; UNIVERSAL; BIO (1-2)

Spekulum Typ:

A; B; C

Spekulum Größe:

S; M; L

Handschuh:

Latex / Nitril

Kollektoren Optionen im Sets:

- Zervikalbürste (I-II)
- Zytologiebürste (auch mit Ball);
- Abstrichspatel Herz-Finger- Form (I-III)
- Bürste und Abstrichspachtel
- Löffel nach Volkmann
- Endorette
- Abstrichtupfer
- Objekträger

Produktnormen:

VM® Standard

Risikoklasse:

Is

Verpackung:

Polybag / Karton

Haltbarkeit:

5 Jahre

EN VM® Gynaecological Kits

- › Are designed for Gynaeco- logical examinations and obstet- ric as well as Gynaecological procedures
- › Are available in various configurations and can be customized according to individual preferences:
 - Vaginal speculum
 - Shoe covers
 - Napkins
 - Gynaecologicalcollectors
 - Examination gloves
 - Condoms
 - Microscope slides

Materials:

PS; PE; PP; Non-wovens; latex/ NBR; nylon; steel; ABS

Sterilization:

Ethylene oxide

Kit type:

MINI; BASIC; CYTO-Standart; CYTO-Premium; ENDO; UNIVERSAL; BIO (1-2)

Speculum type:

A; B; C

Speculum size:

S; M; L

Glove:

Latex / Nitrile

Collector options available in the sets:

- Cervical brush (I-II)
- Cytology brush (also with a ball)
- Ayre spatula (I-III)
- Brush and Spatula
- Volkmann spoon
- Endorette
- Urinary Swab

Compliance standards:

VM® standard

Risk classification:

Is

Packaging:

polybag / carton

Shelf life:

5 years

VM® Gynecological Kits

Index No.	Type	Speculum Size	Gloves Material	Napkin	Microscope slides	Ayre spatula	Urinary swab	Volkman spoon	Endorette	Cytology brush	Cervical brush
-----------	------	---------------	-----------------	--------	-------------------	--------------	--------------	---------------	-----------	----------------	----------------

947	MINI		Latex	•							
-----	------	--	-------	---	--	--	--	--	--	--	--

With speculum type A

948	BASIC	M	Latex	•							
949	CYTO-Standard	M	Latex	•	•	•					
950	CYTO-Premium	M	Latex	•	•	•				•	
951	BIO-1	M	Latex	•			•				
952	BIO-2	S	Latex	•				•			
953	BIO-2	M	Latex	•				•			
954	BIO-2	L	Latex	•				•			
955	BIO-Premium	M	Latex	•	•		•	•			
956	ENDO	M	Latex	•	•				•		•
957	UNIVERSAL	M	Latex	•	•	•	•	•	•	•	•

With speculum type B

958	CYTO-Basic	S	Latex	•		•				•	
959	CYTO-Basic	M	Latex	•		•				•	
960	CYTO-Basic	L	Latex	•		•				•	
961	CYTO-Premium	M	Latex	•	•	•				•	
962	CYTO-Premium	S	Latex	•	•	•				•	
963	CYTO-Premium	L	Latex	•	•	•				•	
964	BASIC	M	Nitrile	•							
965	BASIC	S	Nitrile	•							
966	CYTO-Standard	M	Nitrile	•	•	•					
967	CYTO-Premium	M	Nitrile	•	•	•				•	
968	BIO-1	M	Nitrile	•			•				
969	BIO-2	M	Nitrile	•				•			
970	BIO-Premium	M	Nitrile	•	•		•	•			
971	ENDO	M	Nitrile	•	•				•		•
972	UNIVERSAL	M	Nitrile	•	•	•	•	•	•	•	•

With speculum type C

973	BASIC	M	Latex	•							
974	CYTO-Standard	M	Latex	•	•	•					
975	CYTO-Premium	M	Latex	•	•	•				•	
976	BIO-1	M	Latex	•			•				
977	BIO-2	M	Latex	•				•			
978	BIO-Premium	M	Latex	•	•		•	•			
979	ENDO	M	Latex	•	•				•		•
980	UNIVERSAL	M	Latex	•	•	•	•	•	•	•	•