

PROTECTIVE EQUIPMENT



DE VM® Nitril-Untersuchungshandschuhe

- › Sind vielseitig im klinischen Alltag und schützen das Personal vor gefährlichen und schädlichen Substanzen
- › Sind resistent gegen Thränen, Schnitte, Löcher, Scheuermittel und Chemikalien, welches eine gute Alternative ist, für Personen, die allergisch auf Latex sind

Form:
Beidhändig

Menge Stück im Box / Karton:
100 / 1000; XL: 90/900
200 / 2000; XL: 180/1800

Produktnormen:
EN 455, ASTM D6319,
ISO 11193-1

Materialien:
NBR
Latexfrei

Sterilisation:
Nicht steril

Typ/Länge (mm. min.):
Standard: 240

Externe Oberfläche:
Texturiert

Interne Oberfläche:
Chloriert, Puderfrei

Manschette:
Gleichmäßig gedrehte Walze

Risikoklasse:
I

Verpackung:
Box / Karton

Haltbarkeit:
5 Jahre

EN VM® Nitrile Examination Gloves

- › Are non-sterile medical gloves, intended for conducting medical checkups, diagnostic and therapeutic procedures and for handling contaminated medical material
- › Ideal for medical workers allergic to latex

Shape:
Ambidextrous fitting

Quantity of pieces in box/box:
100 / 1000; XL: 90/900
200 / 2000; XL: 180/1800

Compliance standards:
EN 455, ASTM D6319,
ISO 11193-1

Materials:
NBR
Latexfree

Sterilization:
Non-sterile

Type/length (min. mm):
Standard: 240

External surface:
Textured

Internal surface:
Chlorinated, powder free

Cuff:
Beaded cuff

Risk classification:
I

Packaging:
Box / carton

Shelf life:
5 years



NON-TOXIC. NON-PYROGENIC

Index No.	Size	Glove dimensions (min. mm)		Single wall thickness (min. mm)		Colour	Quantity box/ carton
		Length	Width	Finger	Palm		

Powder free · Finger textured · Non-sterile

981	XS	240	≤ 80	0,07	0,05	Blue	100/1000
982	S	240	80 ± 10	0,07	0,05	Blue	100/1000
983	M	240	95 ± 10	0,07	0,05	Blue	100/1000
984	L	240	110 ± 10	0,07	0,05	Blue	100/1000
985	XL	240	≥ 110	0,07	0,05	Blue	90/900

Powder free · Finger textured · Non-sterile · EKOFREUND (biodegradable)

986	XS	240	≤ 80	0,07	0,05	Pastel Green	100/1000
987	S	240	80 ± 10	0,07	0,05	Pastel Green	100/1000
988	M	240	95 ± 10	0,07	0,05	Pastel Green	100/1000
989	L	240	110 ± 10	0,07	0,05	Pastel Green	100/1000
990	XL	240	≥ 110	0,07	0,05	Pastel Green	90/900

Powder free · Finger textured · Non-sterile

991	XS	240	≤ 80	0,05	0,05	Violet blue	100/1000
992	S	240	80 ± 10	0,05	0,05	Violet blue	100/1000
993	M	240	95 ± 10	0,05	0,05	Violet blue	100/1000
994	L	240	110 ± 10	0,05	0,05	Violet blue	100/1000
995	XL	240	≥ 110	0,05	0,05	Violet blue	90/900

Powder free · Finger textured · Non-sterile

996	XS	240	≤ 80	0,07	0,05	Blue	200/2000
997	S	240	80 ± 10	0,07	0,05	Blue	200/2000
998	M	240	95 ± 10	0,07	0,05	Blue	200/2000
999	L	240	110 ± 10	0,07	0,05	Blue	200/2000
1000	XL	240	≥ 110	0,07	0,05	Blue	180/1800

Powder free · Palm textured · Non-sterile

1001	XS	240	≤ 80	0,10	0,08	Blue	100/1000
1002	S	240	80 ± 10	0,10	0,08	Blue	100/1000
1003	M	240	95 ± 10	0,10	0,08	Blue	100/1000
1004	L	240	110 ± 10	0,10	0,08	Blue	100/000
1005	XL	240	≥ 110	0,10	0,08	Blue	90/900