

# PROTECTIVE EQUIPMENT



## DE VM® OP-Masken

- › Zum Schnüren hinter dem Kopf oder mit einem beidseitigen Elastikband wird aus einem hypoallergenen Vliesmaterial hergestellt und hat hervorragende Atmungsaktivität
- › Bietet einen hohen Widerstand gegen Bakterien, Viren und Staub während der Atmung
- › Eine gute Belüftung verringert mögliche Hautreizungen und ermöglicht ungehindertes Atmen
- › Bietet eine optimale Passform

### Risikoklasse:

I / Is

### Verpackung:

Polybag / Karton

### Materialien:

Vliesstoff

### Sterilisation:

Ethylenoxid

### Typ:

- II
- IIR

### Größe:

Einheitsgröße

### Produktnormen:

EN 13795

## EN VM® Surgical Masks

- › Are made from hypoallergenic non-woven material
- › Offer outstanding breathability and can be secured behind the head, using ties or using an elastic band provided on both sides
- › Ensure a high degree of resistance to bacteria, viruses and dust while breathing
- › Good ventilation reduces potential skin irritation and enables personnel to breathe freely
- › The ergonomically-designed mask has an ideal fit

### Risk classification:

I / Is

### Packaging:

polybag / carton

### Materials:

Non-woven fabric

### Sterilization:

Ethylene oxide

### Type:

- II
- IIR

### Size:

One size

### Compliance standards:

EN 13795

