


# INJECTION MANAGEMENT



## DE VM® Zweiteilige Spritzen

- › Werden zur Verabreichung von flüssigen Medikamenten verwendet
- › Glasklarer Spritzenzylinder mit gut lesbarer Skalierung gewährleistet eine kontrollierte Dosierung
- › Leichtgängiger Spritzenkolben gewährleistet eine schmerzfreie Injektion

**Risikoklasse:**  
ohne Kanüle: Is  
mit Kanüle: IIa

**Verpackung:**  
Blister / Box / Karton

**Haltbarkeit:**  
5 Jahre

**Materialien:**

PP; PE; Stahl; EP-Harze  
Latex und DEHP-frei

**Sterilisation:**

Ethylenoxid

**Volumen (ml):**

2 – 20

**Aufsatz:**

Luer Slip: exzentrisch /  
konzentrisch

**Injektionskanüle:**

optional

**Produktnormen:**

EN ISO 7886, EN ISO 7864

## EN VM® 2-piece Syringes

- › Are suitable for all standard injections by mean of standard techniques
- › The transparent barrel ensures the controlled administration of the medication
- › Smooth-glide plunger ensures painless injection
- › Clearly legible graduation for precise and reliable dosage
- › Secure plunger stop prevents loss of medication

**Risk classification:**  
without needle: Is  
with needle: IIa

**Packaging:**  
blister / box / carton

**Shelf life:**  
5 years

**Materials:**

PP; PE; steel; epoxy resin  
Latex and DEHP-free

**Sterilization:**

Ethylene oxide

**Volume (ml):**

2 – 20

**Taper:**

Luer Slip: eccentric / concentric






**Injection needle:**

optional

**Compliance standards:**

EN ISO 7886, EN ISO 7864

## VM® 2-piece Syringes

Index No.	Volume ml	Taper	Type	Injection needle	Outer diameter		Length mm	Wall	Bevel	Colour code	
					Gauge	mm					
83	2 (2,5)	LS	concentric								
84	5 (6)	LS	eccentric								
85	10 (12)	LS	eccentric								
86	20 (22)	LS	eccentric								
87	2 (2,5)	LS	concentric	side	23G	0,6	30	thin	regular	blue	
88	3	LS	concentric	side	23G	0,6	30	thin	regular	blue	
89	5 (6)	LS	eccentric	side	22G	0,7	40	thin	regular	black	
90	10 (12)	LS	eccentric	side	21G	0,8	40	thin	regular	dark-green	
91	20 (22)	LS	eccentric	side	21G	0,8	40	thin	regular	dark-green	

LL - Luer Lock

LS - Luer Slip