

# INJECTION MANAGEMENT



## DE VM® Großvolumigen Spritzen

- › Für Anwendung in der Urologie, Gastroenterologie oder in der Intensivmedizin
- › Glasklarer Spritzenzylinder mit gut lesbarer Skalierung gewährleistet eine kontrollierte Dosierung
- › Kolbenstopfen mit doppeltem Dichtungsring verhindert das Auslaufen des Medikaments
- › Breites Sortiment an Spritzenansätzen ermöglicht Auswahlmöglichkeit je nach Anwendung

### Materialien:

PP; Gummi; Silikonöl; Stahl;  
EP-Harze  
Latex und DEHP-frei

### Sterilisation:

Ethylenoxid

### Volumen (ml):

50 – 150

### Aufsatz:

- Luer Slip: exzentrisch / konzentrisch
- Luer Lock
- Catheter Tip

### Injektionskanüle:

optional

### Produktnormen:

EN ISO 7886, EN ISO 7864

### Risikoklasse:

ohne Kanüle: Is  
mit Kanüle: IIa

### Verpackung:

Blister / Box / Karton

### Haltbarkeit:

5 Jahre

## EN VM® Large Volume Syringes

- › Are interdisciplinary and applicable in urology, gastroenterology or intensive care
- › The transparent barrel ensures the controlled administration of the medication
- › Smooth-glide plunger ensures painless injection without jerking
- › Secure plunger stop prevents loss of medication

### Materials:

PP; rubber; silicone oil; steel;  
epoxy resin  
Latex and DEHP-free

### Sterilization:

Ethylene oxide

### Volume (ml):

50 – 150

### Taper:

- Luer Slip: eccentric / concentric
- Luer Lock
- Catheter Tip

### Injection needle:

optional

### Compliance standards:

EN ISO 7886, EN ISO 7864

### Risk classification:

without needle: Is  
with needle: IIa







### Packaging:

blister / box / carton

### Shelf life:

5 years

## VM® Large Volume Syringes

Index No.	Volume ml	Taper	Type	Injection needle	Outer diameter Gauge	mm	Length mm	Wall	Bevel	Colour code	
140	50 (60)	LS	eccentric	without							
141	50 (60)	LL	concentric	without							
142	50 (60)	LS	eccentric	side	21G	0,8	40	thin	regular	dark-green	
143	50 (60)	LS	eccentric	side	18G	1,2	40	thin	short	pink	
144	50 (60)	LL	concentric	side	21G	0,8	40	thin	regular	dark-green	
145	50 (60)	LL	concentric	side	18G	1,2	40	thin	short	pink	
146	50 (60)*	CT	concentric	side	18G	1,2	40	thin	short	pink	
147	100*	CT	concentric	side	18G	1,2	40	thin	short	pink	
148	50 (60)	amber LL	concentric	without							
149	50 (60)	CT	concentric	without							
150	100*	CT	concentric	without							
151	120*	CT	concentric	without							
152	150*	CT	concentric	without							
153	200	CT	concentric	without							

\*With Luer adapter

LL - Luer Lock

LS - Luer Slip

CT - Catheter Tip