

> Protective Equipment

PROTECTION COVERALLS

Intended Use: Spray/aerosol-penetration resistant, biohazard-protective coverall for use in EVD patient-isolation units for infection prevention and control against viral penetration.



Type 4:

Chemical protective coverall

Clothing to be worn to protect against spray, airborne solid particulates, radioactive particulates, infective agents with antistatic properties.

Type 5/6:

Coverall to protect against chemical hazards

Clothing to be worn to protect against light spray, liquid aerosol or low pressure, low volume splashes and airborne solid particulates, infective agents, radioactive particulates, clothing with dissipative properties.



PROTECTION COVERALLS

STORAGE

| | |
|-------------|--------------------------|
| Temperature | between -5°C and 40°C |
| Humidity | less than 80% |
| Comply with | Regulation (EU) 2016/425 |

STERILIZATION

| | |
|------------|-------------|
| Method | Non sterile |
| Shelf life | 5 years |



| Catalogue No. | Type | Description | Size (cm) | Material Fabric /gms | Quantity polybag/ carton |
|---------------|-------|--|---|---|---|
| CN4013DT | 4/5/6 | One-piece coverall white colour, with hood, zipper at front opening covered by flap, knitted cuffs, elasticated ankles, hood and waist, heat sealed seams with tape, knitted cuffs: 100% polyester | M: 166-174 L: 174-182 XXL: 182-190 XXXL: 190-198 | Microporous polypropylene + polyethylene film, 65g/m ² white colour. | 1 pcs/package, 50 pcs/ctn. Carton size: 50*32*56cm |
| CN4013DTW | 4/5/6 | Same as CN4013DT with elastic cuffs | | | |
| CN4013DTE | 4/5/6 | Same as CN4013DTW with finger loops | | | |
| CN4013D | 5/6 | One-piece coverall white colour, with hood, zipper at front opening covered by lap, knitted cuffs, elasticated ankles, hood and waist, knitted cuffs: 100% polyester Cut and sewn seams | M: 166-174 L: 174-182 XXL: 182-190 XXXL: 190-198 | Microporous 54% polypropylene + 46% polyethylene, 65g/m ² with white colour. | 1 pcs/package, 25 pcs/ctn. Carton size: 50*30*34cm |
| CN4013E | 5/6 | Same as CN4013D with elastic cuffs | | | |
| CN4013DE | 5/6 | Same as CN4013E with finger loops | | | |
| CN4013DF | 5/6 | Same as CN4013D and CN4013E with shoe cover | | | |
| CN4013DEF | 5/6 | Same as CN4013DE with shoe cover | | | |

| | S | M | L | XL | XXL | XXXL |
|--------|---------|---------|---------|---------|---------|---------|
| height | 158-166 | 166-174 | 174-182 | 182-190 | 190-198 | 198-206 |
| chest | 86-94 | 94-102 | 102-110 | 110-118 | 118-129 | 129-141 |
| waist | 74-82 | 82-90 | 90-98 | 98-106 | 106-117 | 117-129 |



Type 4 standards

Category III

- EN 14605
- EN ISO 13892-1
- EN 13034
- EN 14126
- EN 1073-2
- EN 1149-5

Type 5/6 standards

Category III

- EN 13034
- EN ISO 13892-1
- EN 1149-5
- EN 14126
- EN 1073-2

