

# VM<sup>®</sup> SURGICAL MASKS




**Intended Use:** The surgical masks are made from hypoallergenic non-woven material. It offers a high degree of resistance to bacteria, viruses and dust while breathing.

## Product Description

- Improved filtration
- Without glass fibres
- Hypoallergenic
- Very low resistance to breathing
- Nose bar adaptable
- High filtration capacity
- Perfect fitting

# VM<sup>®</sup> SURGICAL MASKS



## PRODUCT MATERIALS

**Components:** Surgical masks with loops

**Raw Materials:** PP nonwoven 3 ply

## QUALITY AND SAFETY

**Norms** EN 14683

**Shelf life** 5 years

**Risk classification** I

**Efficiency level** Type II

## STORAGE

**Temperature** Between - 15°C and +48°C

**Humidity** less than 80%

## PACKAGING

**Unitary** Box  
Quantity per box: 50

**Carton** Quantity per carton: 2000

## PRODUCT SPECIFICATION

Catalogue No.	Description	Size	Sterilization	Material Fabric	Material/gms	Colour	Quantity box/carton
1314313C	Surgical mask with two soft elastic ear loops	17,5 x9,5 cm	non-sterile	spunbond + meltblown + spunbond	23+25+25	blue	50/2000