

## TECHNICAL DATA SHEET

Ref: CN615

### Description:

CPE Gown, Thumb-loop cuff, Waist ties.

### Intended Use:

Gown designed to protect health professionals and patients from cross infection.

Single-use. Latex free.



### Material Composition:

CPE, 52g ±2g /sqm

**Colour:** Blue

### Packaging:

25pcs/bag, 75 pcs/carton.

Carton size: 38\*20\*31 cm

Gross weight 4.3kg

1\*40'- container – 2760 ctns = 207 000 pcs



### Manufacture:

The quality system of manufacturing sites is ISO13485 compliant.



## TECHNICAL DATA SHEET

### Regulatory information:

Products CE marked as per 93/42/EEC Directive on Medical Devices.

### Storage:

Store in a dry and cool place, away from intense sources of heat.

Keep as much as practicably possible in its shipper box.

### Country of Origin


China

### Shelf life:

3 years shelf life

### Dimension of Finished Product

	XL	Tolerance +/-
A	104.1	3.8
B	111.8	2.5
C	83.8	2.5
D	208.3	5.1
E	31.8	1.3
F	88.9	1.3
G	116.8	2.5
H	68.6	2.5
I	30.5	1.3
J	7.6	1.3
K	15.2	1.3
L	19.1	1.3
M	7.6	1.3
N	7.6	1.3
O	20.3	2.5



*Back View*

